



Email : sales@laboratoryinstrumentindia.com

Product Name :
Air Conditioner (1.5Hp Split Type, Cold & Hot Inverter/
R410a)

Product Code :
LBNY-0002-690003



www.laboratoryinstrumentindia.com

Description :

Air Conditioner (1.5Hp Split Type, Cold & Hot Inverter/ R410a)

Technical Specification :

These heat pump systems are built around Inverter and variable compressor technology, using environmentally friendly R410a refrigerant.

Each mini split indoor unit is designed to heat and cool a single room, and will blend seamlessly into any decor. The unobtrusive front panel features a virtually hidden digital LED display, and the unit operates quietly, day after day.

24 Hours Timer Weekly timer Dry mode & Sleep mode, Multiple speed fan

Low ambient cooling Cold prevention

Compact and Quiet Design I Feel mode, Low voltage start-up.

Air Conditioner (1.5Hp Split Type, Cold & Hot Inverter/ R410a), Air Conditioner (1.5Hp Split Type, Cold & Hot Inverter/ R410a) Equipments, Air Conditioner (1.5Hp Split Type, Cold & Hot Inverter/ R410a) Tools, Air Conditioner (1.5Hp Split Type, Cold & Hot Inverter/ R410a) Tool Kits, Air Conditioner (1.5Hp Split Type, Cold & Hot Inverter/ R410a) Manufacturers, Air Conditioner (1.5Hp Split Type, Cold & Hot Inverter/ R410a) Suppliers from India, China, Kenya

Laboratory Instrument Manufacturers India | Laboratory Instruments Suppliers India

