



Email : sales@laboratoryinstrumentindia.com

Product Name :
Fundamentals of Temperature Measurement

Product Code :
LBNY-0005-103000250



Description :

Experiments in the fundamentals of temperature measurement with 7 typical measuring devices
Various heat sources or storage units: laboratory heater, immersion heater, vacuum flask
Psychrometer for humidity measurement
Tool box for sensors, cables, measuring strips and immersion heater
Mercury, bimetallic and gas pressure thermometers
Calibration units: precision resistors and digital multimeter
Various temperature measuring strips
Temperature sensors: PT100, thermocouple type K, thermistor.

Technical Specification :

Technical Data:

Immersion heater
Power output: 300W
Adjustment of power feed via power-regulated
Socket
Laboratory heater with thermostat
Vacuum flask: 1L.
Power output: 450W
Max. temperature: 425C.



**Laboratory Instrument Manufacturers India | Laboratory Instruments
Suppliers India**