



Email : sales@laboratoryinstrumentindia.com

Product Name :
Air water Heat Transfer Apparatus

Product Code :
LBNY-0005-11200044



www.laboratoryinstrumentindia.com

Description :

The unit is supplied with manuals which describe all parts of the unit, the installation and utilization procedures, as well as many exercises with experimental results. The Air water Heat Transfer Apparatus allows to perform several experiments relevant to the thermal exchange phenomena connected to convection and conduction and some interesting experiments on fluid dynamics. In particular, it is possible to study the thermal exchange between a gas (air) and a liquid (water) and the thermal exchange between an air flow rate in a forced motion and a wall at different temperature.

Technical Specification :

Weight and Dimensions:
Dimensions: 2500x700x1700 (h) mm.
Weight: 190 kg.



**Laboratory Instrument Manufacturers India | Laboratory Instruments
Suppliers India**

