



Email : sales@laboratoryinstrumentindia.com

Product Name :
Four Speed Automatic Transmission

Product Code :
LBNY-0005-12000048



Description :

Electronically controlled.
This transmission designed for high performance has the following parts cut away:
All hydraulic couplings
Torque converter with lockup clutch
Plan.
Electronically controlled. The cut-away parts are:
Planetary gear trains
Differential and hydraulic control
Torque converter with lockup clutch
All hydraulic couplings.
The transmission can easily be turned.

Technical Specification :

Four Speed Automatic Transmission



Laboratory Instrument India

