



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Car Electrical Systems Trainer

**Product Code :**  
LBNY-0005-120000193



**Description :**

Diagnosis of electrical systems is essential at the car service workshop. The basis of diagnosis is the competence of recognising the electrical basic principles applied in the diverse, networked vehicle systems, how they interact in a partial system and how to verify this in measurements. The system comprises all required components; it is compact, can be used at once, and permits free wiring. The focus of the pedagogically valuable course documents is on car diagnosis.

In the end these initially seemingly chaotic partial and overall systems can be reduced to simple circuits where the basic principles apply. When using vehicle lamps in the basic circuits and collecting experience with non-linear behaviour, the learners make their first step into the professional world. The unique course comprises comprehensive documentation and extensive worksheets for teachers and students. The solutions book contains additional further explanations and solution charts. It can be assumed that the trainees targeted here are not dealing with electrical engineering for the first time. Still, the first tasks may commence with basic circuits already known to confirm what was taught before.

**Technical Specification :**

Car Electrical Systems Trainer



**Laboratory Instrument India**

