



Email : sales@laboratoryinstrumentindia.com

Product Name :
Four WD Transfer Case With Locking Differential Cutaway

Product Code :
LBNY-0005-121000325



Description :

It is permanent transfer with possibility of differential locking in order to make the front axle independent from the rear one. Speed reducer transfer case for 4WD cars.

Technical Specification :

Equipment Layout: On Base
Operated: Manually (Through a Crank Handle)
Dimensions: 40 x 40 x 60(h) cm
Net Weight: 25 kg
Gross Weight: 40 kg.



**Laboratory Instrument Manufacturers India | Laboratory Instruments
Suppliers India**