




Email : sales@laboratoryinstrumentindia.com

Product Name :
L Box Stainless

Product Code :
LBNY-0005-12200032



Description :

L Box Stainless steel made consisting of:
Container with inside rigid surfaces,
Two marks far 200 and 400 mm are engraved on the horizontal bottom position.
Obstacle for med by 3 vertical bars having 16 mm dia. with free light of 38 mm,
gate in guillotine form.

Technical Specification :

Dimensions:700 x 300 x 650 mm
Weight :40 kg



Laboratory Instrument India