



Email : sales@laboratoryinstrumentindia.com

Product Name :
Eight Channels Control Unit Touch Screen Semi Automatic

Product Code :
LBNY-0005-12200048



Description :

The new PC based and touch screen control unit developed by department to offer a new generation of control units for the management of semi-automatic and fully automatic compression, flexural, and splitting tensile tests on concrete and cement specimens.

Main Features:

- Modular structure
- Pc based control unit
- No limits memory
- Touch screen display with icons Internet and Intranet connections for remote technical assistance.

Technical Specification :

Eight Channels Control Unit Touch Screen Semi Automatic



Laboratory Instrument India