



Email : sales@laboratoryinstrumentindia.com

Product Name :
Set of 2 Carts

Product Code :
LBNY-0004-780004



Description :

Set of 2 Carts

Technical Specification :

Mechanical bench accessory 02:
Set of 2 carts with:
Built-in force sensors with at least of ± 100 N
3-axis accelerometer
Motion encoder (for measuring position & velocity on or off a track)
Rechargeable battery
Magnetic bumper for force sensor
3-position plunger
Mass tray
Force sensor hook and rubber bumper
Red and blue colors (or any different 2 colors)



Laboratory Instrument India

