



Email : sales@laboratoryinstrumentindia.com

Product Name :
Walnut of Fixing On Bench

Product Code :
LBNY-0004-780010



Description :

Walnut of Fixing On Bench

Technical Specification :

Mechanical bench accessory 08:

Walnut of fixing on bench

Easy Installation: This professional engineering tool can be easily connected to any tool.

The table clamps objects quickly in case of a table for easy grinding, drilling, and carving.

It can hold objects of various regular shapes.

Widely used for carving Bodhi and Walnut, perfect for correcting clocks, sidestands, models, woodworking, and more.



**Laboratory Instrument Manufacturers India | Laboratory Instruments
Suppliers India**

