



Email : sales@laboratoryinstrumentindia.com

Product Name :
Depth Filtration

Product Code :
LBNY-0005-12700035



Description :

Electromagnetic flow rate sensor
4 ball valves with electric drive
Measurement of flow rate, differential pressure, system pressure and temperature
20 tube manometers to measure the pressures in the filter bed
Plotting of Micheau diagrams
Depth filtration with sand filter
Filter velocity adjustable
Sand filter backwash possible
Raw water and backwash pump.

Technical Specification :

Tanks for raw water and treated water
Capacity: each 180L
Raw water pump
Max. head: 10m
Max. flow rate: 13m³/h
Sand filter
Height: 1660mm
Outer diameter: 200mm
Inside diameter: 150mm
Backwash pump
Max. head: 37m

Max. flow rate: 3m³/h.



Laboratory Instrument India