



Email : sales@laboratoryinstrumentindia.com

Product Name :
Fundamentals Of Sedimentation

Product Code :
LBNY-0005-12700039



Description :

5 transparent tanks with scale for comparison of the settling velocities of solids in various suspensions
Experiments in the fundamentals of sedimentation
3 measuring cups for preparation of suspensions
Pycnometer to determine the density of the liquids and solids
Tanks removable for filling, mixing and cleaning
Recommended accessories: balance, coagulant
Tanks backlit by fluorescent tubes to aid observation
Stopwatch to record the sedimentation time.

Technical Specification :

Tanks
Material: PMMA
Length: 1000mm
Scale division: 1mm
Inside diameter: 42mm
Measuring cups
Capacity: 2000ml
Scale division: 50ml
Fluorescent Tubes
Power: 6x 18W
Stopwatch
Resolution: 1/100sec.



**Laboratory Instrument Manufacturers India | Laboratory Instruments
Suppliers India**