



Email : sales@laboratoryinstrumentindia.com

Product Name :
Assembly Exercise Piston Compressor

Product Code :
LBNY-0005-12700051



Description :

Single-stage compressor, air-cooled, with fan flywheel, intake filter and pipe unions
Assembly exercise for engineering training
Compressor consisting of piston and cylinder, housing, driving gear, cylinder cover with valves
Complete, disassembled piston compressor with small parts set and 6 assembly jigs in a storage case.

Technical Specification :

Air-cooled single-cylinder piston compressor
Intake Capacity: 115l/Min
Displacement Volume: 63cm³
Speed: 1850min⁻¹
Cylinder Bore: 50mm
Stroke: 32mm
Max. Pressure: 10bar
Dimensions, assembled, LxWxH: 223x256x314mm
Drive Power Output: 700w.



Laboratory Instrument India

