



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Whirling Of Shafts Apparatus

**Product Code :**  
LBNY-0005-14000035



**Description :**

**Experimental Capability:**

Whirling of a shaft with concentrated masses  
Whirling of a continuous shaft  
Whirling of a long thin shaft supported at its ends.  
Whirling of a shaft with intermediate nodes

**Features:**

Precision electronic tachometer  
Three steel specimen shafts  
Four variable position bearing housings provide nodal points  
Infinitely variable control shaft speed  
Three guard rings to limit whirling deflection.

**Technical Specification :**

Supplied complete with AC Control Unit, tachometer, set of specimen shafts and set of inertia discs and bushes. Comprising a bench mounted light alloy bed 1.8m in length and fitted with high speed DC motor.

---

# Laboratory Instrument India

