



Email : sales@laboratoryinstrumentindia.com

Product Name :
Journal Bearing Demonstration

Product Code :
LBNY-0005-14000046



Description :

To meet propelling demands of our valuable clients, we are presenting ourselves as dependable manufacturers, exporters and suppliers of a wide assortment Journal Bearing Apparatus. Owing to their quality, offered range is highly appreciated around the globe as it is manufactured as per international standards. Assisted by technically advanced facilities, this is designed by our skilled professionals who execute their work in an efficient and precise manner.

Features:

- Precisely designed
- Low maintenance
- Highly durable.

Technical Specification :

- Motor Variable speed Motor
- Control Panel For speed control of motor
- Manometer 16 Tubes
- Bearing Diameter 55 mm
- Journal Diameter 50 mm
- Weights 4 adjustable weights.

Laboratory Instrument India

