



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Pelton Turbine Study Unit Bench Top Model

**Product Code :**  
LBNY-0005-16000017



#### **Description :**

The measurements, which can be performed varying the operation conditions, can be easily compared with data resulting from the theory. The educational Mini-Pelton action turbine allows to study the typical parameters of action turbines. The group is supplied with manuals which describe all parts of the Pelton Turbine Study Unit Bench Top Model unit, the installation and utilisation procedures, as well as many exercises with the relative results. The plant is easy to control and needs a little periodic maintenance.

#### **Technical Specification :**

##### **Required services:**

Electrical supply: 220 V AC single-phase, 50/60 Hz, 1.3 kW

Water supply: tank filling

Weight and dimensions:

Dimensions: 1120x720x850 h mm

Net weight: 95 kg



**Laboratory Instrument India**

