



Email : sales@laboratoryinstrumentindia.com

Product Name :
Automotive Electro Trainer

Product Code :
LBNY-0005-18600045



Description :

Automotive Electro Trainer provides a "hands-on" working knowledge of voltage, current, and resistance with series and parallel circuits. Over 250 of these currently used to train. Automotive Electro Trainer kit allows students to gain a foundation in the principles of electricity and electronics. Graphic symbols are screened into each component's acrylic background. This helps students to understand circuit reading and wiring diagrams. It provides a "hands on" working knowledge of voltage, current, and resistance with series and parallel circuits.

Technical Specification :

Automotive Electro Trainer



**Laboratory Instrument Manufacturers India | Laboratory Instruments
Suppliers India**