



Email : sales@laboratoryinstrumentindia.com

Product Name :
Large Capacity Sample Splitter

Product Code :
LBNY-0005-2020005



Description :

Designed for the reduction of test samples which are too large in volume to be conveniently handled. Complete with two collecting pans. Clam shell hopper: 30 litres capacity. Very sturdily constructed, it is totally Cadmium Plated for rust protection. It handles any material from Sand sizes up to diameter 108 mm. Each Chute bar is 12 mm wide so that openings of 12-24-36-48-60-72-84-96-108 mm are possible.

Technical Specification :

Large Capacity Sample Splitter



Laboratory Instrument India