



Email : sales@laboratoryinstrumentindia.com

Product Name :
Tube, snapcap, conic, 2.0ml, BOX-500

Product Code :
LBNY-0007-22600044



Description :

Intended use:
Various clinical laboratory uses including centrifugation and short-term sample storage.

Technical Specification :

Polypropylene material
2.0ml capacity
Non-sterile
Free of DNA, RNase, DNase & PCR inhibitors
Secure snap cap attached to the tube
Conical base
Graduated
Withstands centrifugation force up to 25,000 x g
Compatible with the standard microcentrifuges.

Weight & volume:
Estimated volume: 6.00dm³
Estimated weight: 0.625kg.

Laboratory Instrument Manufacturers India | Laboratory Instruments Suppliers India

