



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

PCROpticalcaps,0.1/0.2ml,8-strip,BOX-125

**Product Code :**

LBNY-0007-22600052

**Description :**

Cap closure for 0.1ml & 0.2ml PCR tube-strips for use on most standard Thermocyclers, RT PCR platforms, DNA Sequencers.

**Technical Specification :**

8-strip format  
Polypropylene material  
Fits on 0.1ml & 0.2ml tube-strips  
Free of DNA, RNA, DNase, RNase, PCR inhibitors  
Easy to apply and remove  
Non-sterile  
Provides tight seal to minimize evaporation  
Optically clear flat-top profile suitable for real-time PCR.

**Weight & volume:**

Estimated volume: 1.05dm<sup>3</sup>

Estimated weight: 0.028kg.

---

## Laboratory Instrument Manufacturers India | Laboratory Instruments Suppliers India

