



Email : sales@laboratoryinstrumentindia.com

Product Name :
Centrifuge, haematocrit

Product Code :
LBNY-0001-2160001



Description :

Intended for clinical laboratory use; Haematocrit/PCV, Blood separation etc.

Technical Specification :

Automatic brake and lid interlock

Automated speed approx. up to 15000 rpm

Timer 1-99 minutes

Backlit LCD display, digital operable by touch /press allowing responsive operation – even when wearing gloves

Capable of handling capillary tubes.

Power requirements 110-220V/50-60Hz, with voltage surge protection

High sample throughput, multi-safety control, personal program function

Plastic housing, metal lid, deep metal chamber which securely holds the rotor

Rotor holds 24 capillary tubes of 75mm.

Weight & volume:

Estimated volume:0.023m³

Estimated weight: 10kg.

Laboratory Instrument Manufacturers India | Laboratory Instruments Suppliers India

