



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Welding Tong

**Product Code :**  
LBNY-0004-820041



**Description :**

Welding Tong

**Technical Specification :**

Welding Tong

Crafted from the forged steel

provide an ultra-strong grip, ideal balance and long-lasting strength. They can fully articulate and supply exceptional vise gripping strength on the work-piece during heavy hammering.

**High Performance:** Used to easily and properly get the heat treating done, able to hold hot steel firmly on your welding bench, it works well with forging of iron, brass and copper. Perfect for beginner and professional blacksmiths.

**Extremely Versatile:** Specially designed with flat jaws, very versatile and effective in handling wider flat stock sizes. Can better to hold the flat materials while forging, also can hold well on ball bearings and when rolling flat stocks.



**Laboratory Instrument India**

