



Email : sales@laboratoryinstrumentindia.com

Product Name :
Coverall, protection, CatIII, Type3b, M

Product Code :
LBNY-0021-22000107



www.laboratoryinstrumentindia.com

Description :

Disposable liquid-tight biohazard protective coverall, for use in EVD patient isolation units for infection prevention and control against viral penetration.
After use, the coverall is disposed in a biohazard waste container, collected and destroyed. This is also applicable if the coverall is damaged (perforation, etc.).

Technical Specification :

Elasticated hood around face.
Material: Light weight, do not contain rubber/ latex.
Protective seams providing barrier equal to fabric.
Zipper with re-sealable flap protecting leakage through seams.
Elasticated cuffs and ankles.
Sleeves with elasticated thumb loop.
Fabric is Infective agent tested against viral penetration at minimum 1.75kPa
Non-sterile
Antistatic treated on both sides.
Single Use, disposable
Size: M
Each coverall has a stitched-in neck label indicating the type and performance of the suit against the below mentioned standards.
Color: White/ yellow/orange.



Laboratory Instrument India