



Email : sales@laboratoryinstrumentindia.com

Product Name :
Water Basin

Product Code :
LBNY-0004-820049



Description :

Water Basin

Technical Specification :

Water Basin

9-gallon water basin

Screw the stand column in the base tank.

Assemble the washbasin on the top of the stand column.

Insert faucet into transparent tube and fix it on the washbasin, then fix the tube with base tank and stand column.

Fill the liquid soap into the tank and assemble the dispenser.

Screw the flexible drain hose on the washbasin.

Fill fresh water into base tank and step foot pump for handwash.



Laboratory Instrument India