



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Hand Pallet Truck

**Product Code :**  
LBNY-0004-820051



**Description :**

Hand Pallet Truck

**Technical Specification :**

Hand Pallet Truck  
Rubber coated handle  
Epoxy Powder coated body  
Hard chrome plated piston  
Single piece pump with galvanized body  
Large diameter polyurethane steering wheel  
Tandem polyurethane load rollers  
Foot Release / Hand Release  
Entry slide for easy entry in closed pallets.  
High load carrying capacity polyurethane wheels & rollers to protect floors.  
Adjustable Push Rods.  
4mm thick high tensile sheet for extra strenght of forks & body.  
Push rods made from solid square bars instead of pipe.  
Heavy duty additional angle supports on lifting lever frame.  
Working load: 1,500- 1600 kg  
Lift height for fork: (205-250) mm  
Fork height: 85-100 mm  
Fork length mm; 1,000-1200  
Fork outer width mm: 550-600



**Laboratory Instrument Manufacturers India | Laboratory Instruments  
Suppliers India**