



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Portable Photometer, Iron in Foods

**Product Code :**  
LBNY-0022-2200039



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Measuring device equipped with light emitting diode (LED) technology and graphic LCD display. Portable Fluorometer to determine quantity of iron in food. Delivered with accessories. Assembled and packed in a portable case.

**Technical Specification :**

Use: Portable, point of care (POC)  
Light source: LED lamp  
Staff qualification: 1 day training  
Weight: Device: 0.45kg, case incl. accessories and device: 1.9kg



**Laboratory Instrument India**