



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Student Balance 200g

**Product Code :**  
LBNY-0004-10100088



**Description :**

Product is a lever arm double disk balance, there are riding code device, with four accuracy and other weights. The non closed beam is made of aluminum alloy, and the knife is made of steel. Maximum weighing 200g, dividing value 0.02g. The allowable error is + 0.5d (dividing value).

**Technical Specification :**

Student Balance 200g



**Laboratory Instrument Manufacturers India | Laboratory Instruments  
Suppliers India**