



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Transmission, Reflection and Refraction of Light

**Product Code :**  
LBNY-0004-101000307



**Description :**

It is used to do experiments: transmission of light in glass rod, plane mirror reflection of light, verification of refraction law of light, measuring water refractive index, critical angle of refraction of light, reversibility of light in reflection and refraction phenomenon.

**Technical Specification :**

Dimension (L\*W\*H CM): 15\*3\*19



**Laboratory Instrument Manufacturers India | Laboratory Instruments  
Suppliers India**