



Email : sales@laboratoryinstrumentindia.com

Product Name :
Thermomixer Heating Block

Product Code :
LBNY-0001-830014



Description :

Thermomixer Heating Block

Technical Specification :

Ideal for efficient and reliable mixing, heating and cooling established routine applications on tubes and plates.
2D Mix-Control technology: mixes liquids in a controlled, circular movement for fast and reliable mixing in tubes and plates

Anti-spill technology reliably prevents lid wetting and cross-contamination

SmartBlocks provide maximum temperature accuracy and block homogeneity

Temperature settings: 1°C to 100°C

Timer: 15 seconds to 99: 30 h continuous

Maximum speed 3000 rpm

Mixing frequency 300 to 3000rpm

Mixing orbit from 3mm

Different exchangeable thermoblocks (24x0.2ml, 24x0.5ml, 24x1.5ml, 24x2.0ml microplate)

Heating rate maximum 7°C/min (for smartblock plates)

Cooling rate maximum 2.5°C/min between 100°C and room temperature

USB interface

Laboratory Instrument India

