



Email : sales@laboratoryinstrumentindia.com

Product Name :
Capacitance Decade Box

Product Code :
LBNY-0002-620014



Description :

Capacitance Decade Box

Technical Specification :

Decade box with rugged quality and high precision.

Perfect for educational use, R and D, service, production line testing and maintenance.

Unit Features: Sturdy, multi-way binding posts, slide switches add/subtract for desired value, durable ABS plastic housing case.

100 pF to 11uF Capacitance

Capacitance Decade Box, Capacitance Decade Box Equipments, Capacitance Decade Box Tools, Capacitance Decade Box Tool Kits, Capacitance Decade Box Manufacturers, Capacitance Decade Box Suppliers from India, China, Kenya



**Laboratory Instrument Manufacturers India | Laboratory Instruments
Suppliers India**

