



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Non-Sterile Cleanroom Shoe Covers

**Product Code :**  
LBNY-00017-11000105



**Description :**

Seamless, non-skid shoe covers are low-lint for contamination control in environments where wet splash and dry particle protection are needed.

Made with poly-film laminated fabric; double bagged for use in clean manufacturing and quality control environments

Elastic ankle opening for a snug fit

Seamless bottom for better comfort; elastic ankle opening

Low lint for contamination control; Not made with natural rubber latex

Recyclable through the Kimberly-Clark Professional RightCycle Program, which turns used lab coats into raw materials for eco-friendly consumer products

A8 Shoe Covers / Case, 3 Bags of 100 Covers; White

Superior non-skid design has 4x more grip than leading alternatives.

Greener Choice Claims

Environmental benefits include:

Recyclable.

**Technical Specification :**

Non-Sterile Cleanroom Shoe Covers

---

# Laboratory Instrument Manufacturers India | Laboratory Instruments Suppliers India

