



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Sleeves

**Product Code :**  
LBNY-00017-11000190



**Description :**

Material is composed of a breathable polyethylene film/spun bond polypropylene nonwoven laminated fabric. Protective sleeves are constructed of a high quality microporous composite material with a low particle count to meet the requirements of controlled environments.

Static dissipative

Sonic welded seams

Heat dissipating material

Low particle count

Tapered sleeve with elastic at both ends.

**Technical Specification :**

Sleeves



**Laboratory Instrument India**