



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Coveralls, Clean-Processed and Sterile

**Product Code :**  
LBNY-00017-11000316



**Description :**

Made using a patented flash spinning process, provides an inherent barrier to particles, microorganisms and non-hazardous light liquid splash.

Bound seams are covered with garment fabric to reinforce the seam and to reduce the potential for particle penetration

Cleanroom coveralls deliver the best balance of protection, comfort, and durability. Clean-processed and gamma sterilized.

Coverall has elastic openings for tighter fit at wrist and ankle

Raglan sleeve design for greater range of motion and comfort

Bound neck for lower particle shedding

Front zipper closure for easy donning and doffing

White.

**Technical Specification :**

Coverall, Clean-Processed and Sterile



**Laboratory Instrument India**

