



Email : sales@laboratoryinstrumentindia.com

Product Name :
Safe Wear High Performance Lab Coat, Deep Blue

Product Code :
LBNY-00017-11000414



Description :

Made from high-quality 45 g/m2 SMS fabric
Detailed finishing, 3-pocket design and front snap closure
The Safe Wear High Performance Lab Coat, Deep Blue is premium protective apparel made from top-quality materials and designed to provide reliable barrier protection without sacrificing comfort or style.
Each bag contains 10+2 units, for a total of 12 vs. 10 per package from the competition
Available in a wide range of styles, colors and sizes
Comply with recommended OSHA best practices.

Technical Specification :

Safe Wear High Performance Lab Coat, Deep Blue



Laboratory Instrument India