



Email : sales@laboratoryinstrumentindia.com

Product Name :
Light Duty Cut-Resistant Nylon Gloves with Thumbslot

Product Code :
LBNY-00017-11000432



Description :

Light Duty Cut-Resistant Nylon Gloves sleeves are double-layered, made of nylon for soft and comfortable wear experience.

Provides an enhanced protection of wrist

Wrist and elastic band at top provide stability for forearms

Ambidextrous and can be worn by left arm or by right arm

Dirt masking color for prolonged wear

Seamless design

Gauge: 24

For light mechanical protection

Silicone free for use in paint and finish process.

Technical Specification :

Light Duty Cut-Resistant Nylon Gloves with Thumbslot

Suppliers India



Laboratory Instrument Manufacturers India | Laboratory Instruments