



Email : sales@laboratoryinstrumentindia.com

Product Name :
400 Disposable Shirt

Product Code :
LBNY-00017-11000447



Description :

400 provides light-weight inherent barrier protection against hazardous dry particles and aerosols, and nonhazardous light liquid splash.

400 shirts provide the protection and durability you need so you can focus on the task at hand.

Can be used in a variety of applications including general maintenance/operations, asbestos and lead abatement, mold remediation, and environmental cleanups.

Excellent abrasion resistance due its unique, nonwoven material

Technical Specification :

400 Disposable Shirt



Laboratory Instrument India