



Email : sales@laboratoryinstrumentindia.com

Product Name :
Neoprene Sleeves

Product Code :
LBNY-00017-11000479



Description :

This neoprene sleeve offers a broad spectrum of chemical resistance with added protection from snags and punctures.

Excellent resistance to a wide range of chemicals

Economically priced arm protection for practically every application.

Extremely flexible and light weight for level of protection

Durable and cost effective.

Technical Specification :

Neoprene Sleeves



**Laboratory Instrument Manufacturers India | Laboratory Instruments
Suppliers India**