



Email : sales@laboratoryinstrumentindia.com

Product Name :
BCR Cleanroom Socks

Product Code :
LBNY-00017-11000558



Description :

Microdenier technology provides the sock with ideal cool feeling and comfort
Provides microdenier softness for operator comfort
Designed for all-day comfort
Because of the enhanced wicking function, contact discomfort associated with wearing regular cotton or polyester socks is greatly reduced
Reduced odor caused by foot perspiration
Reduced sock degradation over extended wear due to oils and sweat
Long lasting freshness and cooling for operator.

Technical Specification :

BCR Cleanroom Socks



**Laboratory Instrument Manufacturers India | Laboratory Instruments
Suppliers India**

