



Email : sales@laboratoryinstrumentindia.com

Product Name :
Superior Glove Contender Aramid Cut-Resistant Sleeves

Product Code :
LBNY-00017-11000702



Description :

Material: Aramid Fiber
Cuff Style: Knit
For Use With (Application): Metal handling, automotive, Metal fabrication, welding
Construction: Aramid Encased, Silica infused Fiber outer layer with Cool max Inner layer
Size: X-Large
Length (English): 22 in.
Includes: Thumb hole.

Technical Specification :

Superior Glove Contender Aramid Cut-Resistant Sleeves



**Laboratory Instrument Manufacturers India | Laboratory Instruments
Suppliers India**

