



Email : sales@laboratoryinstrumentindia.com

Product Name :
Top Dog 8 mil Die Cut 50 in. Apron

Product Code :
LBNY-00017-11000989



Description :

Ideal for workstations with sharp edges or jobs that require lifting, carrying, or the use of knives or scissors
One-piece construction limits the potential for foreign body contamination and reduces harborage points for pathogens.

Designed to support zoned color-coding programs

Top Dog 8 mil Die Cut 50 in. Apron is reversible and designed to fit comfortably over clothes.

FDA-compliant garments

Thicker 8 mil material significantly increases resistance to punctures and tears

Reusable

Color-coded.

Technical Specification :

Top Dog 8 mil Die Cut 50 in. Apron



Laboratory Instrument India